INTENDED PURPOSE

The first standardised one-stage prothrombin time test was developed by Dr. Armand Quick in 1935. It has now become the basic coagulation screening test for the diagnosis of congenital and acquired deficiencies of clotting factors from the extrinsic pathway (factors IX, V, VII and X)^{1,2}. It is also used for the initial monitoring of oral anticoagulant therapy. This test can be used to assess the protein synthesis capability of the liver in chronic or acute hepatic disorders.

Helena Thromboplastin-Li is of rabbit brain origin and resembles human BCT in its low International Sensitivity Index (ISI). The ISI of Helena Thromboplastin-Li is approximately 1.1 (using manual calibration techniques) and is calibrated against the WHO International Reference Preparation³. Helena Thromboplastin-Li is particularly suited to the monitoring of oral anticoagulant therapy and, in conjunction with the appropriate factor deficient plasma, the measurement of factor activity in the extrinsic pathway. Tissue thromboplastin, in the presence of calcium ions, is an activator which initiates the extrinsic pathway of coagulation. When a mixture of tissue thromboplastin and calcium ions is added to normal citrated plasma, the clotting mechanism is activated, leading to a fibrin clot. If a deficiency exists within the extrinsic pathway, the time required for clot formation will be prolonged depending on the severity of the deficiency.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for *in-vitro* diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION1. **Helena Thromboplastin-Li**

REF 5265 (10 x 2.5ml)

REF 5267 (10 x 5ml)

REF 5269 (10 x 10ml)

Each vial contains rabbit brain thromboplastin suspension. The contents of the vial should be mixed well before and during use. The reagent is a fine suspension of rabbit brain tissue.

2. **Calcium Chloride Solution**

REF 5265 (10 x 2.5ml)

REF 5267 (10 x 5ml)

REF 5269 (10 x 10ml)

Each vial contains 0.025M calcium chloride solution. The reagent is a clear, colourless solution.

3. **Other kit components**

Each kit contains Instructions For Use and lot specific reference values insert.

STORAGE AND SHELF LIFE

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

1. **Helena Thromboplastin-Li**Opened vials are stable 1 month at 2...6°C and 5 days at 15...30°C. The working reagent (mixed with calcium chloride) is stable for 10 days at 2...6°C or 5 days at 15...30°C, AC-4 POS 25 °C = 2hrs. **DO NOT FREEZE**.

Signs of deterioration: Large clumps of particles or changes in expected values may indicate product deterioration.

2. **Calcium-chloride Solution**

Opened vials are stable for 1 month at 2...6°C or 5 days at 15...30°C.

Signs of deterioration: Particulate contamination, cloudiness or a change in expected performance may indicate product deterioration.

ITEMS REQUIRED BUT NOT PROVIDED

REF 5159 Plasma Calibrator Set

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000-3000 x g for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

STEP-BY-STEP PROCEDURE

For accurate INR reporting, it is recommended to determine the laboratory specific ISI of the reagent with the testing system in use. The Helena BioSciences Europe Plasma Calibrator Set (REF 5159) should be used for each new reagent batch. The Helena BioSciences Europe INR Reference Set (REF 5490) should be used to check for shifts in the local ISI which have been noted with changes in laboratory temperature and post instrument servicing, amongst other local variances.

For accurate prothrombin activity reporting (%PT), it is recommended to calibrate the instrument in routine use using the Helena BioSciences Europe PT Calibration Set (REF 5504).

The use of serial dilutions of a reference plasma for the %PT curve is not recommended as this can lead to discrepancies caused by the low fibrinogen in the reference plasma dilutions which are not reflected in patient samples having predominantly normal fibrinogen levels.

1. **Automated Methods, AC-4. (REF 5265 x500 samples, REF 5267 x1000 samples, REF 5269 x2000 samples)**
Refer to the appropriate AC-4 Operators Manual for detailed instructions.
Prepare all reagents as instructed under 'composition'.

START=7s	METHOD=Coag	MIX=0	PATIENT (Pat) vol=50µL Pos=35 (CP)	Def. Plasma (DP) Vol=0µL
INCUB=0s	MATH=Lo gY	CLEAN=0	BUFFER (BUF) Vol=0µL	R1 Vol=0µL
RUNTIME =120s	SENS=2	T-CORR= +45% / -5 sec S-CORR=0	CLEAR (CLR) Vol=0µL	R2 Vol=100µL Pos=25

2. **Manual Methods**
 a) Mix sufficient Helena Thromboplastin-Li and calcium chloride solution to complete the anticipated testing for the day and incubate at 37°C for no more than 8 hours.
 b) Pre-warm 0.1ml of the test plasma at 37°C for 2 minutes.
 c) Add 0.2ml of freshly mixed thromboplastin reagent to the plasma while simultaneously starting a stopwatch.
 d) Note the time for clot formation to the nearest 0.1 seconds.

INTERPRETATION OF RESULTS

Results should be reported to the nearest 0.1 seconds and duplicates should agree within 5% of each other.
INR values can be calculated using the following formula:

$$\text{INR} = (\text{PT Time Patient}) / (\text{Mean Normal PT Time})^{\text{ISI}}$$

%PT values can be interpolated from the calibration graph (%PT of PT Calibration plasmas versus measured clot time), which should be a straight line when plotted on log-log graph paper.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

REF 5186 Norm-Trol 1 10 x 1ml
 REF 5499 Norm-Trol 1 10 x 3ml
 REF 5187 Ab-Trol 2 10 x 1ml
 REF 5183 Ab-Trol 3 10 x 1ml
 REF 5499 INR Reference Set 12 x 0.5ml

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range. This is particularly important for local ISI calibration.

Using the manual (tilt-tube) technique, normal values ranging from 12.5 - 16.5 seconds are typical.

Therapeutic Ranges: The British Society for Haematology⁸ recommend the following INR ranges for various clinical conditions:

- 2.0-2.5 Clinical State Prophylaxis of deep vein thrombosis, including surgery on high risk patients. (2.0-3.0 for hip surgery and fractured femur operations).
 Treatment of deep vein thrombosis. Pulmonary embolism. Systemic embolism. Prevention of venous thrombo-embolism in myocardial infarction. Mitrail stenosis with embolism. Transient ischaemic attacks. Atrial fibrillation.
 3.0-4.5 Recurrent deep vein thrombosis and pulmonary embolism. Arterial disease, including myocardial infarction. Mechanical prosthetic heart valves.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena BioSciences or their representatives as a guideline.

Each laboratory should establish its own performance data.

Reproducibility (using the Behnk RackRotor):

INR	Within assay reproducibility	Clot Time (seconds)	Within Assay CV (%)
n=10	17.0	26.9	1.09
n=10	23.0	2.09	
n=100	23.0	2.38	
n=100	25.6	2.15	
Between assay reproducibility	Clot Time (seconds)	Between Assay CV (%)	
n=3	13.1	2.0	
n=3	27.3	3.3	
n=3	44.0	2.2	
Between batch reproducibility	Clot Time (seconds)	Between Assay CV (%)	
n=3	13.1	2.0	
n=3	27.3	3.3	
n=3	44.0	2.2	

Method Comparison: Comparison of INR values determined using Helena Thromboplastin-Li and two WHO Reference Thromboplastins yielded the following correlation equations:

Thromboplastin-Li = 0.9936 (RTF95) + 0.0129 $r^2 = 0.9834$ n = 79
 Thromboplastin-Li = 0.9589 (BCT441) + 0.1543 $r^2 = 0.9500$ n = 79

BIBLIOGRAPHY

- Quick, A.J. 'A Study of the Coagulation Defect in Hemophilia and Jaundice', Am. J. Med. Sci., 1935; 190: 501.
- Biggs, Rosemary (Ed.), Human Blood Coagulation, Haemostasis and Thrombosis, 2nd Edition, Blackwell Scientific Publications, London, 1976.
- Hirsh, J., Poller, L., Deykin, D., Levine, J. and Dalen, J.E., 'Optimal Therapeutic Range for Oral Anticoagulants', Chest, 1989, 95: 65-115.
- World Health Organisation, Expert Committee on Biological Standards 1984, Technical Series 700/19.
- Poller, L., 'Laboratory Control of Anticoagulant Therapy', Sem. Thromb. Haemostasis, 1986; 12: 13-19.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'The value of plasma calibrants in correcting coagulometer effects on International Normalised Ratios (INR): An international multicentre study', Amer. J. Clin. Pathol., 1995; 103 : 359-365.
- Poller, L., Triplett, D.A., Hirsh, J., Carroll, J. and Clarke, K. 'A comparison of lyophilised artificially depleted plasmas and lyophilised plasmas from warfarin treated patients in correcting for coagulometer effects on International Normalised Ratios (INR)', Amer. J. Clin. Pathol., 1995, 103 : 366-371.
- The British Society for Haematology, 'Guidelines on Oral Anticoagulation: Second Edition', J. Clin. Pathol., 1990; 43 : 177-183.

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