helena | BioSciences Europe

www.helena-biosciences.com

Helena BioSciences Europe
Colima Avenue
Sunderland Enterprise Park
Sunderland
SR5 3XB

tel: +44 (0) 191 549 6064

fax: +44 (0) 191 549 6271

email: info@helena-biosciences.com

Other Helena BioSciences Europe offices:

Helena BioSciences Europe
6 Rue Charles Cros-ZAE
95320 Saint Leu La Foret
France

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INTENDED PURPOSE

For the determination of functional Protein S levels in human plasma by clotting assay.

Protein S is a vitamin K-dependent protein synthesised in the liver, vascular endothelium, and megakaryocytes, which plays an important physiologic role in the Protein C Anticoagulant System^{1,2}. This anticoagulant system is one of the major regulators of haemostasis by inhibiting clot formation and by promoting fibrinolysis. Protein S functions as a cofactor for activated Protein C on the vascular membrane to facilitate the degradation of clotting factors Va and VIIIa, down-regulating clot formation. In normal plasma approximately 40% of Protein S circulates as a free molecule, while 60% is complexed with C4b, a plasma protein of the classical complement pathway³. Only Free Protein S is functionally active and able to bind to activated Protein C, while the complexed form of Protein S is not⁴.

The Helena BioSciences Protein S (clotting) assay determines the amount of functional Protein S. In this assay, dilutions of normal plasma are mixed with protein S depleted plasma. The mixed plasma is then activated by a reagent which contains activated protein C and phospholipid. After activation, clotting is initiated by the addition of calcium chloride. The prolongation of the clotting time is directly related to the concentration of the protein S in the patient plasma.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only - DO NOT INGEST. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg)³ HIV 1 and 2 antibody and HCV antibody; however they should be handled with the same precautions as a human plasma sample.

COMPOSITION

1. Protein S Activator

Each vial contains activated protein C, rabbit brain phospholipid and bulking agents.

Preparation: Reconstitute with 1.0ml purified water. Allow to stand for 20 minutes and mix well before use.

2. Protein S Deficient Plasma

Each vial contains 1.0ml of immunodepleted protein S deficient human plasma.

Preparation: Reconstitute with 1.0ml purified water. Allow to stand for 20 minutes and mix well before use.

3. Other kit components

Each kit contains Instructions For Use.

STORAGE AND SHELF-LIFE

Unopened vials are stable until the given expiry date when stored under conditions indicated on the vial or kit label.

1. Protein S Activator

Reconstituted vials are stable for up to 2 weeks at 2...6°C or 3 hours at 37°C. **DO NOT FREEZE.**

2. Protein S Deficient Plasma

Reconstituted vials are stable for 6 hours at 4°C and 1 month at -20°C if thawed only once.

ITEMS REQUIRED BUT NOT PROVIDED

- Cat. No. 5375 Owrens Veronal Buffer
Cat. No. 5386 25mM calcium chloride solution
Cat. No. 5185 SARP (reference plasma)

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000 - 3000 xg for 15 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

STEP BY STEP PROCEDURE

It is recommended that tests be performed in duplicate.

a) Manual method

1. Reconstitute the required vials of reagent. Place the protein S activator reagent and 25mM calcium chloride solution into a 37°C incubator.
2. Prepare a stock standard by mixing 300µl of reference plasma with 4200µl of Owrens Veronal Buffer.

Use the stock standard immediately to make the following standard dilutions:

Stock	standard volume (µl)	Buffer volume (µl)	Activity %
1	400	0	100
2	300	100	75
3	200	200	50
4	100	300	25
5	0	400	0

3. Prepare patient and control dilutions by adding 50µl of plasma to 700µl of buffer. Keep all standard, patient and control dilutions on ice and test within 30 minutes.
4. To a coagulation cuvette add 100µl protein S deficient plasma and 100µl standard, control or patient dilution.
5. Incubate for 1 minute at 37°C.
6. Add 100µl protein S activator reagent.
7. Incubate for 2 minutes at 37°C.
8. Add 100µl calcium chloride solution and start timer.
9. Note clotting time to the nearest 0.1 seconds.
10. Plot the protein S activity (X-axis) versus clot time (Y-axis) on linear graph paper.
11. Interpolate patient and control values from the straight calibration line.

b) Automated methods

Refer to the Operator's Manual provided by the instrument manufacturer for programming instructions.

INTERPRETATION OF RESULTS

Protein S deficiency, either congenital or acquired, may lead to serious thrombotic events such as thrombophlebitis, deep vein thrombosis, or pulmonary embolism. The prevalence of Protein S deficiency has been estimated to be less than 1 case per 300 in the general population. Two-thirds of patients with a congenital deficiency of Protein S (levels less than 50% of normal) may present with venous thrombosis in young adulthood^{5,6}. In young patients (<35 years) with a history of thrombosis, the prevalence may be as high as 15% to 18%⁷. Acquired Protein S deficiency may be seen during pregnancy, oral contraceptive or oral anticoagulant therapy, liver disease, diabetes mellitus, postoperative complications, septicemia and various inflammatory syndromes⁸. A decreased Protein S activity in plasma may be the result of low concentrations or abnormal function of the Protein S molecule.

A decrease in protein S is associated with increased incidence of thromboembolism¹⁰. A decrease in protein S activity does not necessarily indicate a decrease in concentration. The laboratory diagnosis of Protein S deficiency may require both qualitative (functional) and quantitative (antigen level) determinations¹¹. Helena BioSciences provides both Total and Free Protein S ELISA's (Cat. No. 5286 / 5293).

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.

Helena BioSciences supply the following controls available for use with this product:

Cat. No. 5301 SAC-1

Cat. No. 5302 SAC-2

LIMITATIONS:

Lupus antibody positive samples with prolonged aPTT times may give false high protein S values and should be confirmed by free antigen assay⁹. Borderline or low functional Protein S levels on patients with F-V (Leiden) should be confirmed by free Protein S antigen assay⁹. Patients with high levels (>200%) of protein S activity should be tested at multiple dilutions to obtain accurate results. Corrected protein S levels from at least 2 dilutions must agree. Heparin at levels >2 units/ml may interfere by prolonging clot times and hence increasing apparent protein S values. Patients on warfarin may have false low protein S activity⁹.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

The expected normal range for protein S is gender specific⁹.

Males 85 - 157% Females 65 - 136%

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena BioSciences or their representatives as a guideline. Each laboratory should establish its own performance data.

Within assay reproducibility		
	Protein S Level (%)	Within assay CV (%)
n = 10	94	3.96
n = 8	52	7.21

Between assay reproducibility		
	Protein S Level (%)	Between assay CV (%)
n = 10	82	5.95
n = 8	48	7.53

Accuracy

Comparison of the Helena BioSciences Protein S (Clotting) Kit with another commercially available Protein S method (Free protein S by ELISA) showed well correlated results ($r=0.7781$).

Sensitivity

The method is sensitive to 17.6% Protein S.

Linearity

The method is designed to give a linear standard curve from 17.6%-150% protein S.

BIBLIOGRAPHY

1. DiScipio, R.G. & Davie, E.W., 'Characterisation of Protein S; a gamma-carboxyglutamic acid containing protein from bovine and human plasma.' *Biochemistry*, 1979; 18:899-904.
2. Dahlback, B., 'Purification of human vitamin K-dependent Protein S and its limited proteolysis by thrombin.' *Biochem. J.* 1983; 209:837-846.
3. Walker, F.J., 'Protein S and the Regulation of Activated Protein C.' *Semin. Thromb. Hemost.*, 1984; 10:131-138.
4. Dahlback, B. 'Protein S and C4b-Binding Protein: Components involved in the regulation of the Protein C Anticoagulant System.' *Thromb. Haemost.*, 1991; 66:49-61.
5. Comp, P.C., Nixon, R.R., et al. 'Familial Protein S deficiency is associated with recurrent thrombosis.' *J. Clin. Invest.*, 1984; 74:2082-2088.
6. Walker, .F.J., 'Protein S and Thrombotic Disease.' *PSEBM*, 1992; 200:285-295.
7. Gladson, C.L. et al. 'The frequency of type I heterozygous Protein S and Protein C deficiency in 141 unrelated young patients with venous thrombosis.' *Thromb. Haemost.*, 1988; 59:18-22.
8. Alving, B.M. & Comp, P.C., 'Recent advances in understanding clotting and evaluating patients with recurrent thrombosis.' *Am. J. Obstet. Gynecol.*, 1992; 167:1184-1191.
9. Cooper, P., Cooper, S.M., Smith, J.M. and Makris, M. 'A New Functional Assay of Protein S - An Evaluation' In Press, February 2001.
10. Bertina, R.M., 'Hereditary protein S deficiency' *Haemostasis*, 1985; 15 : 241-245.
11. Comp. P.C. 'Laboratory evaluation of protein S status.' *Semin. Thromb. Haemost.*, 1990; 15 : 177-181.

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BIBLIOGRAPHIE

1. Di Scipio, R. G. & Davie, E. W., 'Characterisation of Protein S; a gamma-carboxyglutamic acid containing protein from bovine and human plasma.' *Biochemistry*, 1979 ; 18 : 899-904.
2. Dahlback, B., 'Purification of human vitamin K-dependent Protein S and its limited proteolysis by thrombin.' *Biochem. J.* 1983 ; 209 : 837-846.
3. Walker, F. J., 'Protein S and the Regulation of Activated Protein C.' *Semin. Thromb. Hemost.*, 1984 ; 10 : 131-138.
4. Dahlback, B. 'Protein S and C4b-Binding Protein: Components involved in the regulation of the Protein C Anticoagulant System.' *Thromb. Haemost.*, 1991 ; 66 : 49-61.
5. Comp, P. C., Nixon, R. R., et al. 'Familial Protein S deficiency is associated with recurrent thrombosis.' *J. Clin. Invest.*, 1984 ; 74 : 2082-2088.
6. Walker, .F. J., 'Protein S and Thrombotic Disease.' *PSEBM*, 1992 ; 200 : 285-295.
7. Gladson, C. L. et al.'The frequency of type I heterozygous Protein S and Protein C deficiency in 141 unrelated young patients with venous thrombosis.' *Thromb. Haemost.*, 1988 ; 59 : 18-22.
8. Alving, B. M. & Comp, P. C., 'Recent advances in understanding clotting and evaluating patients with recurrent thrombosis.' *Am. J. Obstet. Gynecol.*, 1992 ; 167 : 1184-1191.
9. Cooper, P., Cooper, S. M., Smith, J. M. et Makris, M. 'A New Functional Assay of Protein S - An Evaluation' In Press, février 2001.
10. Bertina, R. M., 'Hereditary protein S deficiency' *Haemostasis*, 1985 ; 15 : 241-245.
11. Comp. P. C. 'Laboratory evaluation of protein S status.' *Semin. Thromb. Haemost.*, 1990 ; 15 : 177-181.

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LITERATUR

- DiScipio, R.G. & Davie, E.W., 'Characterisation of Protein S; a gamma-carboxyglutamic acid containing protein from bovine and human plasma.' *Biochemistry*, 1979; 18:899-904.
- Dahlback, B., 'Purification of human vitamin K-dependent Protein S and its limited proteolysis by thrombin.' *Biochem. J.* 1983; 209:837-846.
- Walker, F.J., 'Protein S and the Regulation of Activated Protein C.' *Semin. Thromb. Hemost.*, 1984; 10:131-138.
- Dahlback, B. 'Protein S and C4b-Binding Protein: Components involved in the regulation of the Protein C Anticoagulant System.' *Thromb. Haemost.*, 1991; 66:49-61.
- Comp, P.C., Nixon, R.R., et al. 'Familial Protein S deficiency is associated with recurrent thrombosis.' *J. Clin. Invest.*, 1984; 74:2082-2088.
- Walker, .FJ., 'Protein S and Thrombotic Disease.' *PSEBM*, 1992; 200:285-295.

7. Gladson, C.L. et al. 'The frequency of type I heterozygous Protein S and Protein C deficiency in 141 unrelated young patients with venous thrombosis.' *Thromb. Haemost.*, 1988; 59:18-22.
8. Alving, B.M. & Comp, P.C., 'Recent advances in understanding clotting and evaluating patients with recurrent thrombosis.' *Am. J. Obstet. Gynecol.*, 1992; 167:1184-1191.
9. Cooper, P., Cooper, S.M., Smith, J.M. and Makris, M. 'A New Functional Assay of Protein S - An Evaluation' In Press, February 2001.
10. Bertina, R.M., 'Hereditary protein S deficiency' *Haemostasis*, 1985; 15 : 241-245.
11. Comp. P.C. 'Laboratory evaluation of protein S status.' *Semin. Thromb. Haemost.*, 1990; 15 : 177-181.

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BIBLIOGRAFIA

1. DiScipio, R.G. & Davie, E.W., 'Characterisation of Protein S; a gamma-carboxyglutamic acid containing protein from bovine and human plasma.' *Biochemistry*, 1979; 18:899-904.
2. Dahlback, B., 'Purification of human vitamin K-dependent Protein S and its limited proteolysis by thrombin.' *Biochem. J.* 1983; 209:837-846.
3. Walker, F.J., 'Protein S and the Regulation of Activated Protein C.' *Semin. Thromb. Hemost.*, 1984; 10:131-138.
4. Dahlback, B. 'Protein S and C4b-Binding Protein: Components involved in the regulation of the Protein C Anticoagulant System.', *Thromb. Haemost.*, 1991; 66:49-61.

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5. Comp, P.C., Nixon, R.R., et al. 'Familial Protein S deficiency is associated with recurrent thrombosis.' J. Clin. Invest., 1984; 74:2082-2088.
6. Walker, .F.J., 'Protein S and Thrombotic Disease.' PSEBM, 1992; 200:285-295.
7. Gladson, C.L. et al. 'The frequency of type I heterozygous Protein S and Protein C deficiency in 141 unrelated young patients with venous thrombosis.' Thromb. Haemost., 1988; 59:18-22.
8. Alving, B.M. & Comp, P.C., 'Recent advances in understanding clotting and evaluating patients with recurrent thrombosis.' Am. J. Obstet. Gynecol., 1992; 167:1184-1191.
9. Cooper, P., Cooper, S.M., Smith, J.M. and Makris, M. 'A New Functional Assay of Protein S - An Evaluation' In Press, February 2001.
10. Bertina, R.M., 'Hereditary protein S deficiency' Haemostasis, 1985; 15 : 241-245.
11. Comp, P.C. 'Laboratory evaluation of protein S status.' Semin. Thromb. Haemost., 1990; 15 : 177-181.

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