PCA Ratio Kit

Instructions For Use

REF 5546

Kit Ratio D'APC

Fiche technique

PCA-Ratio-Kit

Anleitung

Kit PCA Ratio

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INTENDED PURPOSE

The Helena BioSciences PCA Ratio Kit is to be used for the detection of APC resistance (factor V Leiden defect).

APC resistance can be measured by the APC ratio of an APTT with the addition of preformed activated protein C divided by the ordinary APTT¹, where ratios above 2.2 indicate normality and ratios below 2.2 indicate the factor V Leiden mutation. PC in the presence of PS may be activated to APC endogenously by PCA contained in a venom fraction from the snake Agkistrodon Contortrix Contortrix. The PCA ratio is similarly obtained by dividing the PCA.APTT by the ordinary APTT². When the patient sample is diluted in factor V depleted plasma, vitamin K dependent factors, intrinsic factors and PC and PS are all normalised. Thus PCA may be used instead of APC, and the PCA ratio becomes a definitive test for APC resistance^{3,4} because APC is generated normally by endogenous PC and PS in the factor V depleted plasma. In the test, plasma is diluted 1/5 in factor V depleted plasma. This is then incubated with an APTT reagent with and without added PCA and the clotting times determined by the addition of CaCl₂. The clotting time ratio of PCA.APTT/ordinary APTT is then determined.

WARNINGS AND PRECAUTIONS

The reagents contained in this kit are for in-vitro diagnostic use only - **DO NOT INGEST**. Wear gloves when handling all kit components. Refer to the product safety data sheet for risk and safety phrases and disposal information. Plasma products have been screened and found negative (unless otherwise stated on the kit box or vial) for the presence of Hepatitis B Antigen (HbsAg) HIV 1 and 2 antibody and HCV antibody; however should be handled with the same precautions as a human plasma sample.

COMPOSITION**1. 25mM Calcium Chloride in Saline**

Each vial contains 8ml of 25mM Calcium Chloride in Saline. The solution is ready for use as packaged.

2. PCA.APTT Reagent

Each vial contains lyophilised PCA.APTT reagent. Reconstitute with 2.0ml purified water. Allow to stand for 20 minutes and mix well before use.

3. APTT Reagent

Each vial contains lyophilised APTT reagent. Reconstitute with 2.0ml purified water. Allow to stand for 20 minutes and mix well before use.

4. Factor V Depleted Plasma

Each vial contains lyophilised Factor V Depleted Plasma. Reconstitute with 2.0ml purified water. Allow to stand for 20 minutes and mix well before use.

5. APC Resistant Plasma

The vial contains lyophilised APC Resistant Plasma. Reconstitute with 0.5ml purified water. Allow to stand for 20 minutes and mix well before use.

6. Other kit components

Each kit contains Instructions For Use.

STORAGE AND SHELF-LIFE

1. 25mM Calcium Chloride in Saline

Calcium Chloride should be stored at 2...6°C, and is stable until the expiry date indicated on the label. AC-4 POS 26 @37°C = 6hrs.

2. PCA.APTT Reagent

Unreconstituted PCA.APTT Reagent should be stored at 2...6°C, and is stable until the expiry date indicated on the label. Reconstituted PCA.APTT Reagent is stable for 2 weeks at 2...6°C. AC-4 POS 32 @15°C = 6hrs.

3. APTT Reagent

Unreconstituted APTT Reagent should be stored at 2...6°C, and is stable until the expiry date indicated on the label. Reconstituted APTT Reagent is stable for 2 weeks at 2...6°C. AC-4 POS 31 @15°C = 6hrs.

4. Factor V Depleted Plasma

Unreconstituted Factor V Depleted Plasma should be stored at 2...6°C, and is stable until the expiry date indicated on the label. Reconstituted Factor V Depleted Plasma should be discarded after use, or can be frozen at -20°C and thawed once.

5. APC Resistant Plasma

Unreconstituted APC Resistant Plasma should be stored at 2...6°C, and is stable until the expiry date indicated on the label. Reconstituted APC Resistant Plasma can be frozen and thawed.

SAMPLE COLLECTION AND PREPARATION

Plastic or siliconised glass should be used throughout. Blood (9 parts) should be collected into 3.2% or 3.8% sodium citrate anticoagulant (1 part). Separate plasma after centrifugation at 2000xg for 10 minutes. Plasma should be kept at 2...6°C. Testing should be completed within 2 hours of sample collection, or plasma can be stored frozen at -20°C for 2 weeks or -70°C for one month. Thaw quickly at 37°C prior to testing. Do not keep at 37°C for more than 5 minutes.

STEP BY STEP PROCEDURE

It is recommended that tests be performed in duplicate.

a) Manual method

1. Prepare duplicate tubes containing 20µl of sample and 80µl of factor V depleted plasma.
2. To one tube add 100µl of APTT reagent and to the duplicate tube 100µl of PCA.APTT reagent.
3. Incubate both tubes for 5 minutes at 37°C.
4. Add 100µl of 25mM CaCl₂/saline to both tubes and record the clotting times.
5. Calculate the PCA.APTT/APTT clotting time ratio.

b) Automated Methods, AC-4. (5546 x80 samples)

Refer to the appropriate AC-4 Operators Manual for detailed instructions.
Prepare all reagents as instructed under 'composition'.

-APC

START=15s	METHOD=Coag	MIX=YES	PATIENT (Pat) vol=10µL pos=36 (CP)	Def. Plasma (DP) Vol=40µL Pos=33
INCUB=300s	MATH=lin	CLEAN=0	BUFFER (BUF) Vol=0µL Pos=0	R1 (APTT) Vol=50µL Pos=31
RUNTIME=300s	SENS=3	T-CORR=0 S-CORR=0	CLEAR (CLR) Vol=0µL	R2 (CaCl ₂) Vol=50µL Pos=26

APCR

START=15s	METHOD=Coag	MIX=YES	PATIENT (Pat) vol=10 μ L pos=36 (CP)	Def. Plasma (DP) Vol=40 μ L Pos=33
INCUB=300s	MATH=lin	CLEAN=0	BUFFER (BUF) Vol=0 μ L Pos=0	R1 (PCA.APTT) Vol=50 μ L Pos=32
RUNTIME=300s	SENS=3	T-CORR=0 S-CORR=0	CLEAR (CLR) Vol=0 μ L	R2 (CaCl ₂) Vol=50 μ L Pos=26

INTERPRETATION OF RESULTS

In a study by the tilt-tube manual method, the PCA ratio kit was shown to have a much greater discrimination between normal plasma samples, factor V Leiden samples and samples from patients on oral anticoagulant therapy than the method employing exogenous APC². Similar results were obtained in the study by Morse and Standen⁵ and in the MDA evaluation study⁶. The method can be used for patients on oral anticoagulant therapy, in patients with high factor VIII levels (during pregnancy), in stored samples with low factor VIII levels and in patients with PC and PS defects. It will also dilute out the effects of lupus anticoagulants and heparin. The main advantage of PCA over preformed APC is the much greater discrimination observed between normal and APC resistant samples (see Performance Characteristics). PCA is also much more stable than APC.

QUALITY CONTROL

Each laboratory should establish a quality control program. Normal and abnormal control plasmas should be tested prior to each batch of patient samples, to ensure satisfactory instrument and operator performance. If controls do not perform as expected, patient results should be considered invalid.


APC Resistant Plasma is supplied as an abnormal control with this kit.

REFERENCE VALUES

Reference values can vary between laboratories depending on the techniques and systems in use. For this reason each laboratory should establish its own normal range.

PERFORMANCE CHARACTERISTICS

The following performance characteristics have been determined by Helena BioSciences or their representatives as a guideline. Each laboratory should establish its own performance data.

BIBLIOGRAPHY

1. Dahlback B, Carlsson M, Svensson PJ. Familial thrombophilia due to a previously unrecognised mechanism characterised by a poor anticoagulant response to activated Protein C. *Proc Natl Sci* 1993;90:1104-08.
2. Denson KWE, Reed SV, Haddon ME, Davidson S, Littlewood TJ. A more discriminating test for APC resistance and a possible screening test to include Protein C and Protein S. *Thromb Res* 1996;81:151-156.
3. Jorquera JL, Montoro JM, Fernandez MA, Aznar J. Modified test for activated Protein C resistance. *Lancet* 1994;334:1162-3.
4. Denson KWE, Reed SV, Haddon ME. The modified APC resistance test. *Thromb Haemost* 1995;74:995.
5. Morse C, Standen G. Specificity of clotting tests for factor V Leiden. *Brit J Haemat* 1996;95:432.
6. Gardiner C, Mackie IJ, Machin SJ, Cooper P, Malia RG, Makris M. 1999. MDA Evaluation Report No 00064 HMSO.

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